

## Goat anti-rabbit IgG-HRP

### 产品信息

---

货 号	ABK0021W	
规 格	100 µg	
储 存 条 件	12 months from manufacturing, storage at -20 °C	
产 品 描 述	<p>Antibody recognizes rabbit immunoglobulin G</p> <p>Polyclonal goat antibody, HRP conjugated. Immunoglobulin G (IgG) is a monomeric immunoglobulin, built of two heavy chains and two light chains. Each IgG has two antigen binding sites.</p>	
溶解缓冲液	10 mM phosphate, pH 7.2, 0.9% NaCl, 50% glycerol with 0.1% ProClin as a preservative	
产 品 应 用	ELISA    Western Blot	
稀 释 比 例	ELISA	1:5,000 - 1:100,000
	Western Blot	1:5,000 - 1:100,000
产 品 浓 度	1.0 mg/ml (+/- 20 %)	
产 品 稳 定 性	TEMPERATURE, TIME	RESULT
	+4 °C, 21 days	OK
	+37 °C, 7 days	OK
Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding.		