

Goat anti-mouse IgG-HRP

产品信息

货号	ABK0020W	
规格	100 µg	
储存条件	-20℃	
产品描述	<p>Antibody recognizes mouse immunoglobulin G</p> <p>Polyclonal goat antibody, HRP conjugated. Immunoglobulin G (IgG) is a monomeric immunoglobulin, built of two heavy chains and two light chains. Each IgG has two antigen binding sites.</p>	
溶解缓冲液	10 mM phosphate, pH 7.2, 0.9% NaCl, 50% glycerol with 0.1% ProClin as a preservative	
产品应用	ELISA Western Blot	
稀释比例	ELISA	1:5,000 - 1:100,000
	Western Blot	1:5,000 - 1:100,000
产品浓度	1.0 mg/ml (+/- 20 %)	
产品稳定性	TEMPERATURE, TIME	RESULT
	+4 °C, 21 days	OK
	+37 °C, 7 days	OK
	Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding.	